



# RemoteOB: A Device to Remotely Control Environmental Conditions and Record Real-Time Behavior of Live Cell

Team member: Ricky Chen, Jo Duan, Arlene Fu, Zoe Yan

Simon Fraser University, Engineering Capstone Project

Date: 09 April 2018

## Contact Info

Company name: HarveStar

Contact person: Ricky Chen

Contact phone number: 604-754-6771

Contact email: sca185@sfu.ca

## Introduction

- Incubator is the necessary equipment for any laboratory conducting cell culture and tissue culture work.
- Advanced live-cell imaging microscopes with onstage chamber have been developed for the purposes of enabling precise control of temperature and gases for time-lapse imaging of live cells under both physiological and non-physiological conditions.

## Problem Definition

Common drawbacks of available competitors (EVOS, ZEISS, WPI, LEICA):

- Not able to remotely control the environmental condition inside the chamber
- Not portable and detachable from microscopes
- Charges at least 70k CAD

## Objectives

RemoteOB serves as an accessory for existing inverted microscopes and it consists of four Separable modules:

### Environmental control unit

Detect and adjust current environmental conditions such as *temperature*, *humidity* and *CO2 concentration* in the chamber.

### Environmental onstage chamber

The enclosed chamber where live cell culture is placed.

### Position Adjustment System

Adjust the position between the objective of the microscope and the chamber.

### Imaging System

Record the image or video captured by a digital camera associated with the inverting microscope and transmit the information to web-interface.





Figure 1: System Block Diagram

## Web Interface

### Welcome to RemoteOB



Help Sign Up Log in



### Current Environmental Setting

Temperature

0°C

Humidity

5%

CO2 Concentration

0%

### Target Environmental Setting

Temperature

38

Humidity

5

CO2 Concentration

X-Y Position

Left

Forward

Backward

Right

## Schematic Detail




Figure 2: Blueprint of RemoteOB




Figure 3: Onstage Chamber




Figure 4: X-Y Stage Schematics




Figure 5: Chamber Schematics

## Conclusion

- By finishing the stage of the proof-of-concept, a confident conclusion can be made that the final product can be delivered with promised features
- By assembling RemoteOB to the existing inverted microscopes, the whole integrated system enables remote control of environmental condition and real-time imaging of live-cell growth.

## Results

- Temperature sensor should function properly within designed circuits
- The temperature should be adjusted by 12V heating pad controlled by 5V relay
- The circuit design should consider components for safety and protection purposes
- The dummy values of carbon dioxide and humidity are fed into the system, illustrating the working mechanism of the system feedback
- Data reading from the temperature sensor to the microcontroller
- Enable the two-way data transmission between microcontroller, Raspberry Pi and Web server.
- A basic web interface displays outputs from the temperature sensor and CO2 sensor

## Future Work

- Be able to inject CO2 and N2 into the chamber through a long tube by using the electric valves
- Build the x-y stage position adjustment system
- Build the humidity system by injecting mist into the chamber through the tube that connects to Ultrasonic mist-maker to provide enough humidity for the cell to grow.
- Integrate existing modules and planned modules together and test the overall functionality of control unit and web app

## Acknowledge

Daniel Cao, Product Development Engineer, Evo-in-motion Technology, EIM PCB

Contact email: zcao1991@gmail.com

## Reference

ThermoFisher, "EVOS Onstage Incubator," [Online]. Available: <https://www.thermoFisher.com/ca/en/home/life-science/cell-analysis/cellular-imaging/cell-imaging-systems/evos-l-auto/evos-onstage-incubator.html>. [Accessed 26 March 2018].

Leica, "Inverted Microscope Solution DMI8 S Platform," Leica, [Online]. Available: <https://www.leicamicrosystems.com/products/light-microscopes/inverted-microscopes/details/product/leica-dmi8/>. [Accessed 25 March 2018].